

# BREEZI

A Range of Practical  
Paediatric Chairs



TEEZI



BREEZI



BREEZI MAX

# BREEZI RANGE



Built and assembled in the UK, Breezi range of chairs are neat pieces of craftsmanship that can be incorporated into any home, school, college or nursery.

Every chair can be used as a simple seating solution, or made bespoke to meet the needs of each child.



They are easily adjustable, meaning they can grow with the child



A wide range of seat depths and optional extras are available to create the optimal seating position



Various mounted castor options for easy manoeuvrability at home or school



The Breezi range caters for all sizes from tiny tots to young adults



Upholstered accessories are available in a variety of colours and fabrics, with the option to personalise



Eco-friendly wooden design is easy to clean, and can be incorporated amongst home furnishings



The Breezi Range can be used to help children with mild to moderate needs, including those with postural and behavioural difficulties.

## Measurements of the Breezi family

	Width	Depth	Height
Teezi	43 cm	54 cm	68 cm
Breezi	46 cm	62 cm	78 cm
Breezi Max	51 cm	69 cm	93 cm

## Head Positioning



Simple Flat Head Rest



Head Rest with Wings

## Upper Body Positioning



Chest Strap

## Pelvic Positioning




Hip Belt



Thigh Pads



Knee Abductor (Pommel)



Pre-ischial Ridge Cushion

## Lower Body Positioning



Heel Bar



Foot Shoes



Corner Back

## Back Options



Standard Back Rails



High Back Rails



High Back Rails (inverted)



High Back Seat

### Mobile Options



Millie-Mova



Push Handles



Skis with Castors



Low Level Mobile Frame



Mobile Base



Trays



Stabilising Skis



Arm Rest



Grab Rails

### Upholstery

Cushioning is available for a range of upholstered items. Here is a list of all items available with cushioning:

**Seat Cushion**

**Standard Back Rail Pad**

**High Back Rail Pad**

**High Back Pad**

**Cornerback Padding**

**Thigh Pads**

**Pommel**

**Tray Pads**

**Head Rest Pad**

**Pre-ischial Ridge Cushion**

**Chest Strap**

**Pelvic Trap**